
630
ART. 114.630

GranberG®

630
ART. 114.630

GranberG®

Available sizes: XS-5/6, S-6/7, M-7/8, L-8/9, XL-9/10
Tilgjengelige størrelser: XS-5/6, S-6/7, M-7/8, L-8/9, XL-9/10
Tillgängliga storlekar: XS-5/6, S-6/7, M-7/8, L-8/9, XL-9/10
Saatavilla olevat koot: XS-5/6, S-6/7, M-7/8, L-8/9, XL-9/10
Dostępne rozmiary: XS-5/6, S-6/7, M-7/8, L-8/9, XL-9/10

EN: USER MANUAL

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SV: BRUKSANVISNING

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GranberG®

Head office:
GRANBERG AS
Bjoavegen 1442, NO-5584 Bjoa, Norway
Tel +47 53 77 53 00, e-mail post@granberg.no

Swedish office / European warehouse:
GRANBERG SVERIGE AB
Schubergvägen 20, 311 74 Falkenberg, Sweden
Tel +46 (0)346 124 25, e-mail post@granberg-ab.se

granberggloves.com

ENGLISH

Single-Use Gloves Magic Touch®. Soft Nitrile™, powder-free. Violet blue colour.

Medical Device Regulation

- This product is classified as a Class I medical device and meets the provisions of the Medical Device Regulation (EU) 2017/745.
- This product complies with European Standards EN 455-1:2020, EN 455-2:2015, EN 455-3:2015, and EN 455-4:2009.

EU Type-Examination

- This product is classed as Category III Personal Protective Equipment (PPE) according to European PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 and EN ISO 374-5:2016.
- Notified Body responsible for certification and Module B compliance is SATRA Technology Europe Ltd., Bracetown Business Park, Clonee, Dublin D15 YN2P, Ireland. Notified Body number 2777.
- Notified Body for internal production control plus supervised product checks at random intervals (Module C2) is SATRA Technology Europe Ltd., Bracetown Business Park, Clonee, Dublin D15 YN2P, Ireland. Notified Body number 2777.

EN ISO 374-1:2016 permeation levels are based on breakthrough times as follows:

Permeation performance level	1	2	3	4	5	6
Measured breakthrough time (min)	>10	>30	>60	>120	>240	>480

Chemical Hazards Pictogram: EN 16523-1:2015; Additional information on chemical resistance obtainable from manufacturer. This product complies with type C requirements and the following pictogram shall be used with reference to clause 6.3 of ISO 374-1.

Micro Organism Hazards Pictogram: EN ISO 374-5:2016 - Protection against Bacteria, Fungi and Virus. If there is no penetration of bacteriophages through the specimen, the following pictogram is applied.

AQL is 1.5; Additional information obtainable from the manufacturer.

Gloves in compliance with the Regulation (EC) 1935:2004 and according to Commission Regulation (EU) No. 10/2011 suitable for handling all types of foods for a short period of time.

Performance and Limitation of Use

- This product has been tested in accordance with EN ISO 374-5:2016. Protection against bacteria and fungus - Pass. Protection against viruses - Pass.
- Gloves had been tested in accordance with EN 16523-1:2015 resistance to permeation by chemicals and in accordance with EN 374-4:2013:

	EN ISO 374-1:2016 (Type C)	Level	EN 374-4:2013 Degradation, %
K	40% Sodium Hydroxide	6	-9.5
T	37% Formaldehyde	4	16.1

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm - where the cuff is tested also) and relates only to the chemical tested. The result may differ if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defect or imperfections.

EN 374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemicals.

- This product provides protection against Bacteria, Fungi and Virus. The gloves had been tested in accordance with ISO 16604:2004 (the test method is uses Phi-X174 bacteriophages) to meet the requirements of BS EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
- Components used in glove manufacturing may cause allergic reactions in some users. If allergic reactions occur, seek for medical advice immediately.
- These gloves are shorter than the minimum length requirements stated in EN 420 as they are designed for special purpose where minimum risk at wrist area and for examination purpose.

Product Instruction for Use

- Usage – For Single Use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes; and increased risk of holes and tear during re-use due to weakening of gloves by cleaning processes.
- Sizing – Select the right size glove for your hand.
- Donning – Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to a get a good fit. Don the other glove by the same procedure.
- Inspection – Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.
- During Use - Gloves polluted by blood or other organic fluids must be immediately removed as medical waste. In case of pricking or tearing, user should immediately change gloves.
- In case of feeling looseness on fingers, user should change gloves.
- Doffing – Hold glove bead and pull toward the finger until the glove come off.

Storage Conditions

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight.

Ideally stored in dry conditions between 10 °C to 30 °C in its original package. Do not exceed 40 °C.

When stored as recommended, the gloves will not suffer any changes in chemical properties for up to three years from date of manufacture. Service life cannot be specified and depends on application and responsibility of user to ascertain suitability of the glove for its intended use.

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