



Product description

Order number	85.1017
Product description	Safety lancet, Extra, Ø needle: 18 G, yellow, penetration depth: 1.8 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

Product characteristics

Application	capillary blood collection
Example of use	Children, Adults, obese adults: capillary blood gas analysis, CBC
Puncture site	finger, ear, heel
Penetration depth	1.8 mm
Safety product	yes

Size

Blood volume	20 µl - 100 µl
Diameter of needle	18 G
Width of product	21 mm
Length of product	52 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Safety lancet, Extra, Ø needle: 18 G, penetration
depth: 1.8 mm**

Page 2

Material & colours

Colour of product yellow

Purity & certification

Product category Medical device
CE certified CE 0124 (DEKRA)
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1200
Type of inner packaging case
Piece(s) / inner box 200
Piece(s) / case 1200
Piece(s) / pallet 108000
Depth of box 118 mm
Width of box 118 mm
Height of box 151 mm
Depth of case 366 mm
Width of case 248 mm
Height of case 173 mm
Case volume 0.0157 cbm
Weight of product 0.0029 kg
Weight of case 4.08 kg
EAN of inner box 4038917085045
EAN case 4038917057646

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-06-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT